gibthai
A 3N HOLDING COMPANY

MOLECULAR DIAGNOSIS

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
TB-Diagnosis

WHO ORGANIZATION

2.1 Conventional tests for the diagnosis of TB

In many high TB burden settings, sputum smear microscopy remains the primary diagnostic technique for evaluating individuals who present with the signs and symptoms of TB, or screen positive for TB. However, sputum smear microscopy is an insensitive test, with a limit of detection (LoD) of 5000–10 000 bacilli/mL of sputum (16). Furthermore, sputum smear microscopy cannot distinguish drug-susceptible strains from drug-resistant strains. WHO continues to recommend that NTPs replace microscopy as the initial diagnostic test with WRDs that directly detect MTBC.

WRD
NTP

WHO-recommended rapid diagnostic test
national TB programme

WHO operational handbook on tuberculosis

Module 3: Diagnosis

2025 Update

TB-LAMP

WHO ORGANIZATION

World Health
Organization

WHO's policy recommendations

1. TB-LAMP may be used as a replacement test for sputum-smear microscopy to diagnose pulmonary TB in adults with signs and symptoms consistent with TB (conditional recommendation, very low-quality evidence).
2. TB-LAMP may be used as a follow-on test to smear microscopy in adults with signs and symptoms consistent with pulmonary TB, especially when further testing of sputum smear-negative specimens is necessary (conditional recommendation, very low-quality evidence).

WHO consolidated guidelines on tuberculosis

Module 3: Diagnosis

Rapid diagnostics for
tuberculosis detection

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
Mycobacterium Tuberculosis

Nucleic Acid Test card

One-Step Molecular Testing Solution

Pluslife Mini Dock

PM001

Rapid Molecular POCT System

Pluslife Dock Pro 8

PM008

Eight-Channel POC Molecular Analyzer

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


Loop mediated isothermal amplification

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
Independent Developed Isothermal amplification Technique: RHAM

RHAM (Isothermal)

The Solution of lacking Sensitivity and Specificity

- The RHAM system includes a probe reporter system consisting of specially modified probes and nucleases.
- A unique probe reporter system further amplifies the signal to ensure all signals are derived from specific amplification products, making RHAM superior in sensitivity and specificity to LAMP.



Analytical Performance

Sensitivity (Limit of Detection)

The analytical sensitivity of the assay was performed using three lots of reagents tested across three testing days with known number of *M. tuberculosis* H37Rv. The limit of detection is the lowest concentration of organisms at which 19 of 20 replicates were positive. **The analytical sensitivity is 50 CFU/mL.**

The test cover many species of MTB including

- *Mycobacterium tuberculosis*
- *Mycobacterium bovis*
- *Mycobacterium africanum*
- *Mycobacterium canetti*
- *Mycobacterium caprae*
- *Mycobacterium pinnipedii*

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Sputum 0.5 ml

Mixing: 1 – 2 mins hand on time
with 10 mins sample incubation

75°C

Result: 10 ~ 30 mins

Sputum

Three Products with WHO Endorsement

	Cepheid	Molbio	Pluslife
Sample	Sputum	Sputum & Tongue swab	Sputum & Tongue swab
Pre-processing	Liquefaction treatment	External nucleic acid extraction	Extraction-free
Duration	~ 120 min	~ 60 min	5 - 25 min
Device			

Pluslife is much faster, cheaper and more convenient than Cepheid & Molbio, while delivering similar performance

[Health Topics ▾](#)[Countries ▾](#)[Newsroom ▾](#)[Emergencies ▾](#)[Data ▾](#)[About WHO ▾](#)

Department for HIV, Tuberculosis, Hepatitis and Sexually Transmitted Infections

[< Back to Near point-of-care tests, tongue swabs, and sputum pooling for TB](#)

NPOC-NAATs

Near point-of-care nucleic acid amplification tests (NPOC- NAATs) as a new diagnostic class for diagnosis of TB using sputum and tongue swabs

Near point-of-care nucleic acid amplification tests (NPOC-NAATs) are swab-based molecular tests for TB detection that can produce results from a primary sputum or tongue swab sample in less than one hour using instruments that can be battery operated and do not require specialized infrastructure for use or storage.

These tests can be done in basic peripheral laboratories, such as those that perform smear microscopy, and health clinics, mobile units, or community sites that do not have laboratories. They can be performed by health care workers with basic technical skills because they do not require laboratory methods like precision pipetting. However, like smear microscopy, this class of tests does not incorporate rifampicin-resistance detection and therefore requires reflex testing with a complementary solution for drug-resistance determination.

More information on the defining characteristics of TB NPOC-NAATs may be found [here](#).

Recommendation

In adults and adolescents with signs and symptoms of pulmonary TB or who screen positive for pulmonary TB, NPOC-NAATs on sputum should be used as initial diagnostic tests for TB rather than smear microscopy (*strong recommendation, moderate certainty of evidence*);

When sputum cannot be obtained, NPOC-NAATs on tongue swabs should be used as initial diagnostic tests for TB (*strong recommendation, moderate certainty of evidence*).

Clinical Study Report of MTBC Nucleic Acid Test Card (Sputum)

MTBC Test Result Summary (n=99)

Pluslife	Xpert MTB/RIF Ultra		
	Positive	Negative	Total
Positive	71	0	71
Negative	1	27	28
Total	72	27	99

Study Purpose

The purpose of the study is to evaluate the performance of the MTBC Nucleic Acid Test Card for detection of Mycobacterium Tuberculosis complex DNA in the **raw sputum** samples as compared to results determined by Xpert MTB/RIF ultra as a reference test.

Kit performance

- Sensitivity = 98.61 % (95% CI: 92.54 to 99.75%)
- Specificity = 100 % (95% CI: 87.54 to 100%)
- Accuracy = 98.99 % (95% CI: 94.50 to 99.82%)

■ Kappa score = 0.9748

Conclusion

The study results demonstrated **good clinical performance** of the MTBC Nucleic Acid Test Card to detect Mycobacterium tuberculosis complex (MTBC) DNA in raw sputum samples from patients who are clinically suspected of pulmonary tuberculosis.

PQ: MiniDock MTB Test Matches GeneXpert® Ultra, Exceeds WHO TPP, and Outperforms Smear Microscopy

Diagnostic accuracy against MRS, comparator tests (N=1266*)			
Test (Sample)	Sensitivity n/N (% , 95% CI)	Specificity n/N (% , 95% CI)	p-value [†]
MiniDock MTB Test (Sputum Swab)	189/218 (86.7% , 81.5–90.9)	1026/1048 (98.0% , 97.0–98.8)	0.11
Xpert Ultra (Sputum)	195/218 (89.4% , 84.6–93.2)	1038/1048 (99.0% , 98.3–99.5)	
Test (Sample)	Sensitivity n/N (% , 95% CI)	Specificity n/N (% , 95% CI)	p-value [†]
Smear Microscopy (Sputum)	136/218 (62.4% , 55.6–68.8)	1046/1048 (99.8% , 99.3–100)	<0.001

* Includes only those with valid results on all four index tests and valid MRS results; Excludes participants with sputum Xpert Trace results

^ Includes 95 participants with a single negative smear

† For sensitivity difference based on McNemar's test of paired proportions

Reference: R2D2 TB Network. Clinical Performance Study of MiniDock® MTB Test in High TB Burden Settings. Stage Analysis Report, 14 Jul 2025.

MiniDock MTB Tests (Sputum Swab) **matches** GeneXpert® Ultra in **both sensitivity and specificity** and **far better than** sputum smear microscopy.

Diagnostic Performance of the Pluslife MiniDock MTB and Molbio MTB Ultima Assays to Detect Tuberculosis From Tongue and Sputum Swabs Among Outpatients and in Active Case Finding in Cameroon

Cyrille Mbuli,^{1,2} Jean Bosco Taya-Fokou,^{1,2} Rita Nsamenang,^{1,2} Nestor Biatu,^{1,2} Gildas Nguimfack,^{1,2} Ngha Ndze Mbuh,^{1,2} Joceline Konso,^{1,2} Zourriyah Adamou Mana,^{1,2} Arthur Rausmis Nankouo-Njouonang,^{1,2} Magdalene Tiamuh,^{1,2} Denis Nsame,² Norah Nyah Ndi,^{3,4} Irene Adeline Goupeyou Wandji,^{4,5} Mercy Fundoh,⁵ Maurice Ganava,^{6,7} Ousmanou Bello,⁷ Paul Meoto,^{8,9} Adeline Fitame,⁹ Comfort Vuchas,^{1,2} Teyim Pride Mbuh,^{10,11} Valerie Flore Donkeng-Donfack,^{11,12} Tushar Garg,¹² Jacob Creswell,¹² and Melissa Sander^{1,2}; the RAPID TB Consortium^a

¹Center for Health Promotion and Research, Bamenda, Cameroon; ²Bamenda Regional Hospital, Bamenda, Cameroon; ³Cameroon Baptist Convention Health Services, Bamenda, Cameroon; ⁴National TB Program—Littoral Region, Douala, Cameroon; ⁵National TB Program—Northwest Region, Bamenda, Cameroon; ⁶National TB Program—Far North Region, Maroua, Cameroon; ⁷National TB Program—North Region, Garoua, Cameroon; ⁸National TB Program—Southwest Region, Buea, Cameroon; ⁹National TB Program—West Region, Bafoussam, Cameroon; ¹⁰Tuberculosis Reference Laboratory Douala, Douala, Cameroon; ¹¹Centre Pasteur du Cameroun, Yaoundé, Center, Cameroon; and ¹²Stop TB Partnership, Geneva, Switzerland

Diagnostic Accuracy comparator tests (N= 1097)		
Test	Sensitivity (%, 95% CI)	
	culture-positive TB	smear-negative TB
MiniDock MTB Test	86% (95% confidence interval [CI], 79.0–91.0)	67% (95% CI, 51.0–79.0)
Molbio MTB Ultima	84% (95% CI, 76–89)	58% (95% CI, 41.0–74.0)

Sputum*

Sputum: Pluslife vs Cepheid ~ 100%

1

Source: Italy

Culture or Cepheid			
Pluslife	Positive	Negative	Total
Positive	13	0	13
Negative	0	2	2
Total	13	2	15

Sensitivity: 100% Specificity: 100% Agreement: 100%

2

Source: Kenya

Culture or Cepheid			
Pluslife	Positive	Negative	Total
Positive	38	0	38
Negative	0	24	24
Total	38	24	62

Sensitivity: 100% Specificity: 100% Agreement: 100%

3

Source: India

GOVERNMENT OF KARNATAKA
HEALTH & FAMILY WELFARE SERVICE
OFFICE OF THE JOINT DIRECTOR(TB), LADY WILKINGDON, STATE TB CENTRE,
2nd FLOOR, ARDRA KUDURU, MAGADI ROAD, BENGALURU-560023

No.LWSTC/NTEP/STDC/101/2023-24 Date:04.10.2023

Performance evaluation of Huel UniAMP Mycobacterium Tuberculosis Nucleic Acid Test Card

Product name: Mycobacterium Tuberculosis Nucleic Acid Test Card
Trade name: UniAMP MTB Card

Description
Mycobacterium tuberculosis is a pathogenic bacterial species belonging to the genus Mycobacterium, and is the causative agent of Tuberculosis. Tuberculosis has been a leading cause of death, especially for low- and middle-income countries. The rapid and accurate diagnosis of Mycobacterium tuberculosis complex (MTB-complex) infection is important to reduce morbidity and mortality rates and risk of transmission.

Intended use
The Mycobacterium Tuberculosis Nucleic Acid Test Card is used for in vitro qualitative detection of Mycobacterium Tuberculosis complex DNA in human pulmonary and extrapulmonary samples. It is recommended to combine the patient's clinical manifestation and other laboratory tests to analyse the condition.

Testing principle
The assay is based on isothermal amplification method and enzyme digestion probe technology, and conserved regions are selected for specific primers and specific probe design. A large number of target sequence's copies were generated in the reaction system during the isothermal amplification. When the probe hybridizes to the complementary sequence, it is cleaved and fluorescence is emitted. Integrated Nucleic Acid Testing Device detects and analyses fluorescence signal automatically, reporting negative, positive or invalid result. The assay includes internal control for monitoring of sample collection, processing, and amplification to reduce false negative results.

Objective:
1. Ease of operation in comparison to TrueNat and GeneXpertMTB Rif assay.
2. Comparative performance evaluation of UniAMP MTB (Huel) against GeneXpertMTB Rif assay(Cepheid).

Methodology:
Sputum Swab:
① With clean hands, slightly extract the swab without touching its delicate tip.
② Carefully swirl the swab 10 times by dipping the bristles inside the sputum/EP sample except tissue sample collected in the container for 10 -15 seconds

DEPARTMENT OF CLINICAL MICROBIOLOGY
CHRISTIAN MEDICAL COLLEGE
Ida Scudder Road, Vellore - 632 004, India.

Dr. V. Balaji, MD, PhD, FRC Path, Professor
Dr. Jay S Michael, MD, FRC Path, Professor
Dr. Asha A John Prekshith, MD, Professor
Dr. Ravi Dhanu Sekar, MD, DTM&S Professor
Dr. Bhaskar Lal S, MD, Professor

Dr. Ravi Subramanian P, MD, FRC Path, Associate Professor
Dr. Manjree M Nair, DCP, MD, Associate Professor
Dr. Sanketh Karanavara P, MD, DTM&S, Associate Professor
Dr. S. Dharmadheeran, MD, Assistant Professor
Dr. Lydie Jeyaseelan S, MD, Associate Professor

Performance evaluation of Huel UniAMP Mycobacterium Tuberculosis Nucleic Acid Test Card

Name of the kit: UniAMP Mycobacterium tuberculosis nucleic acid test card
Batch number: PQLMTB2301, PQLMTB2302, PQLMTB2303
Manufacturing date: Aug 2023
Expiry date: Aug 2024

Please find below the evaluation report of the following Mycobacterium tuberculosis nucleic acid kit.

	CB NAAT		Total
	Detected	Not detected	
UniAMP MTB	Detected	34	38
	Not Detected	2	42
	Total	36	77

Conclusion:
1. Based on the methodology
• UniAMP test procedure does not require decontamination or extraction of DNA.
• Test procedure takes 7-15 min in comparison to CB-NAAT and TrueNat.
• Very easy to operate doesn't require dedicated facility with specific low room temperature
• UniAMP can be implemented in field, a true POCT device.
2. Based on test performance
CB-NAAT TB detection kit and UniAMP MTB kit performances are comparable. Out of 77 samples tested 38 are positive and 42 samples were negative by UniAMP MTB. Three samples were discordant, of these 2 samples which were negative and a single sample which was positive could not be repeated due to inadequate sample.

Dr. Jay S Michael, MD, FRC Path, Professor
Department of Clinical Microbiology
Christian Medical College,
Vellore - 632 004, T.N., India

Sputum	CB NAAT		
Pluslife	Detected	Not Detected	Total
Detected	33	0	33
Not Detected	0	7	7
Total	33	7	40

Sensitivity: 100% Specificity: 100% Agreement: 100%

Source: GOVERNMENT OF KARNATAKA,HEALTH &FAMILY WELFARE SERVICE


Sputum	CB NAAT		
Pluslife	Detected	Not Detected	Total
Detected	27	0	27
Not Detected	1	32	33
Total	28	32	60

Sensitivity: 96.4% Specificity: 100% Agreement: 98.3%

Source: DEPARTMENT OF CLINICAL MICROBIOLOGY


Pluslife MTB Solution Backed by Leading Partners: Multi-Center Trials in 10+ Countries

Studies

- A multi-center clinical study across **>10** LMICs
- **>5,000** Samples, sponsored by NGOs
- Collaborating with NGOs such as R2D2, GHL, and FIND

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Division of Tuberculosis validation

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No.	Sample LSN.	Xpert result	Pluslife result	Result
1	RAM 67-8770	Not detected	Negative	Pass
2	WHH 67-8775	Not detected	Negative	Pass
3	PAN 67-8782	Not detected	Negative	Pass
4	WHH 67-8776	MTB Detected (very low)	Positive	Pass
5	LAN 67-8887	MTB Detected (low)	Positive	Pass
6	MJA 67-8939	MTB Detected (low)	Positive	Pass
7	TBC 67-8774	MTB Detected (medium)	Positive	Pass
8	MJA 67-8938	MTB Detected (medium)	Positive	Pass
9	FAH 67-8300	MTB Detected (high)	Positive	Pass
10	SKP 67-8227	MTB Detected (high)	Positive	Pass

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No.	Sample LSN.	Pluslife result			Result
		Test 1	Test 2	Test 3	
1	RAM 67-8770	Negative	Negative	Negative	Pass
2	WHH 67-8775	Negative	Negative	Negative	Pass
3	PAN 67-8782	Negative	Negative	Negative	Pass
4	FAH 67-8300	Positive	Positive	Positive	Pass
5	TBC 67-8774	Positive	Positive	Positive	Pass
6	SKP 67-8227	Positive	Positive	Positive	Pass

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